

Cloning - Promega pGEM-T Easy kit

Ligation

1. Add to the reaction mixture:
 - a. 3 μ l Fresh PCR product
 - b. 5 μ l 2x Rapid Ligation Buffer
 - c. 1 μ l pGEM-T vector
 - d. 1 μ l T4 DNA Ligase
2. Incubate at room temperature for 1 hour or overnight at 4°C
3. Store on ice until needed

Transformation

1. Get cells out of freezer, thaw on ice for 5 minutes
2. Add 2 μ l of ligation mixture into a 1.5 ml tube on ice
3. **Carefully** transfer 50 μ l of cells into the tube
4. Incubate on ice for 20 minutes
5. Heat-shock the cells for 45-50 seconds at 42°C
6. Immediately put the cells on ice for 2 minutes to recover
7. Add 950 μ l of SOC recovery media
8. Shake at 37°C for 1.5 hours
9. Spread 100 μ l onto LB plates with X-gal and ampicillin
10. Incubate plates overnight at 37°C and you should get colonies.

Materials required

Cleaned PCR product
Promega pGEM-T easy cloning kit
Promega JM109 Competent cells
SOC media
LB + amp + IPTG + X-gal plates
10 μ l, 200 μ l, and 1000 μ l pipette tips
1.5 ml microcentrifuge tubes
Ice

Equipment necessary

10 μ l micropipette

200 μ l micropipette

1000 μ l micropipette

37° C incubator

37° C shaking incubator

42° C heat block