

Cloning – Invitrogen TOPO-TA kit

Ligation

1. Add to the reaction mixture:
 - 4 μ l Fresh PCR product
 - 1 μ l Salt Solution
 - 1 μ l TOPO vector
2. Incubate at room temperature for 5 minutes
3. Store on ice until needed, or overnight at 4° C for maximum number of transformants

Transformation

1. Get cells out of freezer, thaw on ice for 5 minutes
2. Add 2 μ l of ligation mixture
3. Incubate on ice for 5 minutes
4. Heat-shock the cells for 30 seconds at 42° C
5. Immediately put the cells on ice
6. Add 250 μ l of SOC recovery media
7. Shake at 37° C for 1 hour
8. Spread 50 μ l on LB plates with X-gal and ampicillin
9. Incubate plates overnight at 37° C and you should get colonies
10. Proceed to colony PCR to verify your results

Materials necessary

Cleaned PCR product
TOPO-TA cloning system
TOP-10 competent cells
SOC media
LB + amp + IPTG + X-gal plates
10 μ l pipette tips
1000 μ l pipette tips
1.5 ml microcentrifuge tube
Ice

Equipment necessary

37° C shaking incubator
37° C incubator
10 μ l micropipette
1000 μ l micropipette
42° C heat block