

# **Qiagen DNeasy DNA extraction protocol for bacterial cultures**

*Adapted from QIAgen DNeasy handbook, July, 2006.*

## **Procedure:**

1. Appropriately label a 1.5 ml tube for each sample.
2. Add 1.75 ml of bacterial culture to a labeled 2 ml tube.
3. Spin tubes at 20,000 x g for 5 minutes in centrifuge. Decant liquid.
4. Add 180 ul of enzymatic lysis buffer to you tube and vortex 10-20 s.
5. Incubate at 37° C for 30 min.
6. Add 25 ul of proteinase K to the tube
7. Add 200 ul of Buffer AL to the tube.
8. Vortex the tube briefly.
9. Incubate at 56° C for 30 min. *Now is a good time to label all the tubes you need for the rest of the protocol.*
10. Add 200 ul of 100% ethanol to the tube
11. Vortex briefly.
12. Using a micropipette, transfer entire contents (~600 ul) of tube to labeled spin column.
13. Centrifuge column at 10,000 x g for 1 min.
14. Remove column from collection tube. Place column in new collection tube.
15. Add 500 ul of buffer AW1 to the column and centrifuge at 10,000 x g for 1 minute.
16. Remove column from collection tube. Place column in new collection tube.
17. Add 500 ul of buffer AW2 to the column and centrifuge at 20,000 x g for 3 minute.

18. Carefully remove tubes from centrifuge, do not let flow-through contact column. If this happens, spin tube again for 1 min at 20,000 x g.

19. Transfer the column to a 1.5 ml tube and add 200 ul of buffer AE to the column. Let column stand at room temperature for 1 minute.

20. Centrifuge at 10,000 x g for 1 minute. Discard the column and store the DNA appropriately (4° C for short term, -20° C for long term).

**Materials required**

Qiagen DNeasy Blood and Tissue kit  
200 and 1000 ul pipette tips  
1.5 ml microcentrifuge tubes  
2.0 ml microcentrifuge tubes  
Overnight bacterial cultures

**Equipment required**

Bench top centrifuge capable of 20,000 x g  
200 ul micropipette  
1000 ul micropipette  
vortexer