

Squish Buffer DNA extraction

Procedure

1. Homogenize tissue in 50 μ l of squish buffer.
2. Add an additional 100 μ l of squish buffer.
3. Add 1 μ l of Proteinase K (4 μ g/ml).
4. Incubate 60 min at 37° C; increase temperature to 85° C and incubate for 10 min.
5. Use 1-2 μ l of DNA extraction for PCR.

Materials necessary

- 1.5 ml microcentrifuge pestle tubes
- Mini pestles for 1.5ml tubes
- Squish buffer
- Proteinase K
- 10, 20 and 200 μ l pipette tips

Equipment necessary

- 37° C incubator
- 85° C incubator
- 10 μ l micropipette
- 20 μ l micropipette
- 200 μ l micropipette

Note: Squish buffer extractions do not last as long as other methods and should be stored at -20° C for no more than 30 days.

Adapted from the method of Gloor et al. 1991. Targeted gene replacement in *Drosophila* via P element-induced gap repair. *Science* 253: 1110.